

Appendix 21. Cold chain equipment technical sheets

N°	REFRIGERATOR
Brand:	Model:
Net storage capacity:	litres
Gross storage capacity:	litres
Holdover time without power at 43°C:	hours
Drugs in stock:	
Presence of temperature monitoring tools:	
Date:	Thermometer <input type="checkbox"/> Temperature monitoring sheet <input type="checkbox"/> Stop!Watch® <input type="checkbox"/>

N°	FREEZER
Brand:	Model:
Gross storage capacity:	litres = icepacks
Freezing capacity:	kg/24h = icepacks
Holdover time without power at 43°C:	hours
Icepacks in stock:	0.6 litre <input type="checkbox"/> 0.4 litre <input type="checkbox"/> Other <input type="checkbox"/>
Presence of temperature monitoring tools:	
Date:	Thermometer <input type="checkbox"/> Temperature monitoring sheet <input type="checkbox"/>

N°	FRIDGE/FREEZER
Brand:	Model:
Fridge net storage capacity:	litres
Freezer net storage capacity:	litres
Freezing capacity:	kg/24h = icepacks
Holdover time without power at 43°C:	hours
Drugs in stock:	
Icepacks in stock:	0.6 litre <input type="checkbox"/> 0.4 litre <input type="checkbox"/> Other <input type="checkbox"/>
Presence of temperature monitoring tools:	
Date:	Thermometer <input type="checkbox"/> Temperature monitoring sheet <input type="checkbox"/> Stop!Watch® <input type="checkbox"/>